


For in-house culture

**Bringing antimicrobial stewardship into VET practice.
Use the right antibiotic the first time.**




NG-TEST® Vet ESBL

- Rapid detection of CTX-M groups 1, 2, 8, 9, 25 β -lactamases
- From pure bacterial colonies
- User friendly
- Long shelf-life

For use in cats and dogs

NGB-CTM-S23-100 | Kit of 20 tests




NG-TEST® Vet CRE

- Rapid detection and identification of the 5 major types of carbapenemase enzymes (NDM, IMP, VIM, OXA-48 and KPC) in Enterobacterales or aeruginosa
- From pure bacterial colonies
- User friendly
- Long shelf-life


For use in cats and dogs

NGB-CAR-S23-100 | Kit of 20 tests



In the exam room

Save time in sample preparation to optimize the blood culture workflow.



NG-TEST® Vet Direct ESBL

- **Ready-to-use kit** comprising a sample preparation device (DirecTool) and a rapid test in cassette format (CTX-M Multi) designed for the detection of the main extended-spectrum β -lactamases (ESBLs)
- **Simple preparation protocol**
- **Concentration of bacteria** from clinical samples
- **Save valuable time** - No need to wait for culture
- **Adapted for urine and rectal swab**

For use in cats and dogs

NGB-CTM-U23-100

Kit of 10 tests



NG-TEST® Vet Direct CRE

- **Ready-to-use kit** comprising a sample preparation device (DirecTool) and a rapid test in cassette format (CARBA-5) designed for the detection and differentiation of the five most common carbapenemase families
- **Simple preparation protocol**
- **Concentration of bacteria** from clinical samples
- **Save valuable time** - No need to wait for culture
- **Adapted for urine and rectal swab**

For use in cats and dogs

NGB-CAR-U23-100

Kit of 10 tests