


Detection and differentiation of resistance mechanisms

For professional use only

Fast identification of AMR-related enzymes to optimize patient management.



NG-TEST® CARBA-5®


- **Rapid detection and differentiation** of the 5 most common Carbapenemases families (NDM, IMP, VIM, OXA-48-like and KPC) including their most clinically relevant variants produced by Enterobacterales or *Pseudomonas aeruginosa*
- **From bacterial samples** - Validated with NG-TEST® Blood Culture PREP
- **Excellent correlation with PCR**
- **User friendly & long shelf-life**

Sensitivity: 100% - Specificity: 100%

NGB-CAR-S23-021

Kit of 20 tests

CE FDA cleared



NG-TEST® CTX-M Multi

- **Rapid detection** of the 5 most prevalent ESBL CTX-M Groups (1-2-8-9-25) produced by Enterobacterales
- **From bacterial samples** - Validated with NG-TEST® Blood Culture PREP
- **Excellent correlation with PCR**
- **User friendly & long shelf-life**

Sensitivity: 100% - Specificity: 100%

NGB-CTM-S23-016

Kit of 20 tests

CE FDA cleared



NG-TEST® MCR-1

- **Rapid detection** of MCR-1 involved in Colistin resistance
- **From bacterial colony**
- **Excellent correlation with PCR**
- **User friendly & long shelf-life**

Sensitivity: 100% - Specificity: 98.4%

NGB-MCR-S23-016

Kit of 20 tests

CE

Detection and differentiation of resistance mechanisms

For professional use only




NG-TEST® *Acineto-5*

- **Rapid detection and differentiation** of OXA-24-like/143-like, OXA-23-like, OXA-58-like, VIM and NDM carbapenemases from *Acinetobacter spp.* bacterial sample
- **From bacterial samples**
- **Excellent correlation with PCR**
- **User friendly & long shelf-life**

NGB-ACI-S23-000

Kit of 20 tests

RUO



Sample preparation kits

For professional use only

Streamline sample preparation to save time and optimize the blood culture workflow.



NG-TEST® Blood Culture PREP

- Easy and fast preparation of bacterial sample from a positive blood culture
- Validated in combination with NG-TEST® CARBA-5® and NG-TEST® CTX-M Multi
- Enables faster detection and differentiation of ESBL and major carbapenemases

NGP-BCP-S73-001

Kit of 20 tests



NG-TEST® DirecTool

- Concentration of bacteria from samples
- Saves valuable time - No need to wait for culture
- Maximizes laboratory workflow

NGA-DTL-U93-000

Kit of 10 tests